The backpack that will grow if needed. Perfect for a 2–3 day mission







#### **FEATURES**

- The backpack's size is adjustable, ranging from 30 to 40 L.
- The shoulder pads include a chest strap and attachment for a handset.
- There is a carrying handle at the top and on one side of the bag.
- The backpack is equipped with PALS/MOLLE webbing on the outside of the pockets and on the lower sides.
- The main compartment's zipper pullers are made of webbing, while the outside pockets feature a plastic cord end, allowing them to be easily distinguished in the dark.
- Inside the main compartment, there are two pockets, and the interior is lined with Velcro loop fabric, enabling the attachment of additional pouches.
- The outer pockets are zippered and contain mesh pockets inside.
- The backpack supports the use of the spoon system.
- At the top, a zipper can be opened for access or utilized as a port for antennas, cables, or a hydration system hose.

## **VERSIONS**

29-01410A01-000 Black 29-01410A09-000 Grey

#### **DIMENSIONS**

H=50, W=34 & D=17-24 cm. 30-40L

### WEIGHT

1500 g.

# **MATERIAL**

Fabric: 100% polyamide PU coated Webbing: 100% polyester. The fabric is IRR and FR.







- \* Opening for antennae, cables or a hydration tube.
- \* The main compartment's four zipper pullers are made of webbing, while the outside pockets feature a plastic cord end, allowing them to be easily distinguished in the dark.
- \* The exterior pockets are zipped and they have a zipped mesh pocket and three compartment elastic mesh pocket on the inside. MOLLE bands on the outside and two attachment loops for a repair buckle.
- \* Side handle
- \* The main zipper i S-shaped to make it easy to just open a part of the backpack to access it.
- MOLLE webbing on the lower sides so additional pouches can be attached on the backpack.
- All straps have Velcro roll up ends to secure them.



- \* Top handle
- Attachment loops for a repair buckle or side tension strap
- \* Shoulder pads with holder for a microphone.
- \* Chest strap
- \* Padded and ventilated back.
- \* Quick release buckles on each shoulder pad.
- \* Holes for the Spoon sleeve/Curved long Spoon.
- \* Hip strap with side release buckle.



On the top there is a zipped opening for a radio antenna or a hydration system.



The zipped opening in the top can also be used to insert a bolt cutter.



The exterior pockets each have one interior zipped compartment in mesh and one three slot compartment in elastic mesh.

There are two loops in the lower pocket, blue to be easy to find, to secure gear.




If needed, the backpack can be expanded from 30L to 40L capacity. Just open the zipper along the back.



The frame is strong and shaped to match the curve of your back.



The back is shaped after your back to sit comfortably and to make it easier to carry body armour.





The back and front of the inside is covered with Velcro loop fabric. Inside pouches can be attached (not included).



Snigels light combat belt can be used as a padded waist belt. Just insert it behind the mesh/fabric in the lower part of the back. The size can be adjusted by placing it in different heights.



You can easily open just a part of the main zipper and access the backpack where ever you want.



Snigels Spoon system can be used in this backpack. It gives you the possibility to use the hip belt as a combat belt and to ditch the backpack fast. Spoon sleeve, Curved long sleeve and a Spoon garage or double spoon garage.



The Spoon sleeve with Curved long Spoon is inserted in the back of the backpack, in one of the slots. Adjust it to the proper height.



The Spoon sleeve and Curved long Spoon attached for a short/regular length back.



The Spoon sleeve and Curved long Spoon attached for a long length back



The Spoon garage/double Spoon garage is attached on the middle of the hip belt. It can also be used on any MOLLE combat belt.



The Curved long Spoon is inserted into the Spoon garage.



The hip belt can move with your hips and when you ditch the backpack, the hip belt with pouches stays on your body.